### Dr. Jesús Pato Moruelo

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1. **Branemark PI, Hansson BO, Adell R, et al.** Osseointegrated implants in the treatment of the edentulous jaw. Experience from a 10-year period. *Scand J Plast Reconstr Surg Suppl.* 1977; 16: 1-132.
2. **Schulte W, Kleineikenscheidt H, Lindner K, Schareyka R.** The Tübingen immediate implant in clinical studies. *Dtsch Zahnärztl Z.* 1978; 33: 348-359.
3. **Del Fabbro M, Boggian C, Taschieri S.** Immediate implant placement into fresh extraction sites with chronic periapical pathologic features combined with plasma rich in growth factors: Preliminary results of single-cohort study. *J Oral Maxillofac Surg* 2009; 67: 2476-2484.
4. **Quayle AA, Cawood JI, Smith GA, Eldridge DJ, Howell RA.** The immediate or delayed replacement of teeth by permucosal intra-osseous implants: The Tübingen implant system. Part 2: Surgical and restorative techniques. *Br Dent J* 1989; 166: 403-410.
5. **Hammerle CH, Chen ST, Wilson TG Jr.** Consensus statements and recommended clinical procedures regarding the placement of implants in extraction sockets. *Int J Oral Maxillofac Implants* 2004; 19 Suppl: 26-28.
6. **Esposito MA, Koukouloupoulou A, Coulthard P, Worthington HV.** Interventions for replacing missing teeth: Dental implants in fresh extraction sockets (immediate, immediate-delayed and delayed implants). *Cochrane Database Syst Rev.* 2006; (4): CD005968.
7. **Attieh MA, Al Sabeeha NH, Duncan WJ, et al.** Immediate single implant restorations in mandibular molar extraction sockets: A controlled clinical trial. *Clin Oral Implants Res* 2013; 24: 484-496.
8. **Saadoun AP, Missika P, Dennis L.** Immediate placement of an implant after extraction: Indications and surgical requirements. *Actual Odontostomatol (Paris).* 1990; 44: 415-435.
9. **Missika P.** Immediate placement of an implant after extraction. *Int J Dent Symp.* 1994; 2: 42-45.
10. **Annibaldi S, Bignozzi I, Iacovazzi L, La Monaca G, Cristali MP.** Immediate, early and late implant placement in first molar sites: A retrospective case series. *Int J Oral Maxillofac Implants* 2011; 26: 1108-1122.
11. **Cornelliani R, Cangini F, Covani U, Wilson TG Jr.** Immediate restoration of implants placed into fresh extraction sockets for single-tooth replacement: A prospective clinical study. *Int J Periodontics Restorative Dent* 2005; 25: 439-447.
12. **Bersani E, Coppede AR, de Paula Pinto Prata HH.** Immediate loading of implants placed in fresh extraction sockets in the molar area with flapless and graftless procedures: A case series. *Int J Periodontics Restorative Dent* 2010; 30: 291-299.
13. **El-Chaar ES.** Immediate placement and provisionalization of implant-supported, single-tooth restorations: A retrospective study. *Int J Periodontics Restorative Dent* 2011; 31: 409-419.
14. **Becker W, Goldstein M.** Immediate implant placement: Treatment planning and surgical steps for successful outcome. *Periodontol* 2000 2008; 47: 79-89.
15. **Blus C, Szmukler-Moncler S.** Atraumatic tooth extraction and immediate implant placement with piezosurgery: Evaluation of 40 sites after at least 1 year of loading. *Int J Periodontics Restorative Dent.* 2010; 30: 355-363.
16. **Araujo MG, Wennstrom JL, Lindhe J.** Modeling of the buccal and lingual bone walls of fresh extraction sites following implant installation. *Clin Oral Implants Res* 2006; 17: 606-614.
17. **Moya-Villaescusa MJ, Sanchez-Perez A.** Measurement of ridge alterations following tooth removal: A radiographic study in humans. *Clin Oral Implants Res* 2010; 21: 237-242.
18. **Tan WL, Wong TL, Wong MC, Lang NP.** A systematic review of post-extraction alveolar hard and soft tissue dimensional changes in humans. *Clin Oral Implants Res* 2012; 23 (Suppl 5): 1-21.
19. **Chen ST, Wilson TG, Jr, Hammerle CH.** Immediate or early placement of implants following tooth extraction: Review of biologic basis, clinical procedures, and outcomes. *Int J Oral Maxillofac Implants* 2004; 19 Suppl: 12-25.
20. **Becker W, Sennerby L, Bedrossian E, Becker BE, Lucchini JP.** Implant stability measurements for implants placed at the time of extraction: A cohort, prospective clinical trial. *J Periodontol* 2005; 76: 391-397.
21. **Koh RU, Rudek I, Wang HL.** Immediate implant placement: Positives and negatives. *Implant Dent* 2010; 19: 98-108.
22. **Hammerle CH, Araujo MG, Simion M, Osteology Consensus Group 2011.** Evidence-based knowledge on the biology and treatment of extraction sockets. *Clin Oral Implants Res.* 2012;23 Suppl 5: 80-82.
23. **Kois JC.** Predictable single-tooth peri-implant esthetics: Five diagnostic keys. *Compend Contin Educ Dent* 2004; 25: 895-6, 898, 900 passim; quiz 906-7.
24. **Waasdorp JA, Evian CI, Mandracchia M.** Immediate placement of implants into infected sites: A systematic review of the literature. *J Periodontol* 2010; 81: 801-808.
25. **Schropp L, Isidor F.** Timing of implant placement relative to tooth extraction. *J Oral Rehabil* 2008; 35 (Suppl 1): 33-43.
26. **Novaes AB Jr, Novaes AB.** Immediate implants placed into infected sites: A clinical report. *Int J Oral Maxillofac Implants* 1995; 10: 609-613.
27. **Dos Santos MV, Elias CN, Cavalcanti Lima JH.** The effects of superficial roughness and design on the primary stability of dental implants. *Clin Implant Dent Relat Res* 2011; 13: 215-223.
28. **Quiryren M, Van Assche N, Botticelli D, Berglundh T.** How does the timing of implant placement to extraction affect outcome? *Int J Oral Maxillofac Implants* 2007; 22 (Suppl): 203-223.
29. **Vandamme K, Naert I, Geris L, Sloten JV, Puers R, Duyck J.** Histodynamics of bone tissue formation around immediately loaded cylindrical implants in the rabbit. *Clin Oral Implants Res* 2007; 18: 471-480.
30. **Akkocaoglu M, Uysal S, Tekdemir I, Akca K, Cehreli MC.** Implant design and intraosseous stability of immediately placed implants: A human cadaver study. *Clin Oral Implants Res* 2005; 16: 202-209.
31. **Neugebauer J, Weinlander M, Lekovic V, von Berg KH, Zoeller JE.** Mechanical stability of immediately loaded implants with various surfaces and designs: A pilot study in dogs. *Int J Oral Maxillofac Implants* 2009; 24: 1083-1092.
32. **Novaes AB Jr, Papalexou V, Grisi MF, Souza SS, Taba M Jr, Kajiwara JK.** Influence of implant microstructure on the osseointegration of immediate implants placed in periodontally infected sites. A histomorphometric study in dogs. *Clin Oral Implants Res* 2004; 15: 34-43.
33. **Ferguson SJ, Langhoff JD, Voelter K, et al.** Biomechanical comparison of different surface modifications for dental implants. *Int J Oral Maxillofac Implants* 2008; 23:1037-1046.
34. **Grandi T, Garuti G, Samarani R, Guazzi P, Forabosco A.** Immediate loading of single post-extractive implants in the anterior maxilla: 12-month results from a multicenter clinical study. *J Oral Implantol* 2012; 38: 477-484.
35. **Crespi R, Cappare P, Gherlone E.** A 4- year evaluation of the peri-implant placement parameters of immediately loaded implants placed in fresh extraction sockets. *J Periodontol* 2010; 81: 378-83.
36. **Covani U, Orlando B, D'Ambrosio A, Sabattini VB, Barone A.** Immediate rehabilitation of completely edentulous jaws with fixed prostheses supported by implants placed into fresh extraction sockets and in healed sites: A 4-year clinical evaluation. *Implant Dent* 2012; 21: 272-279.
37. **Ganeles J, Zollner A, Jackowski J, ten Bruggenkate C, Beagle J, Guerra F.** Immediate and early loading of Straumann implants with a chemically modified surface (SLActive) in the posterior mandible and maxilla: 1-year results from a prospective multicenter study. *Clin Oral Implants Res* 2008; 19: 1119-1128.
38. **Groisman M, Frossard WM, Ferrelira HM, de Menezes Filho LM, Touati B.** Single-tooth implants in the maxillary incisor region with immediate provisionalization: 2-year prospective study. *Pract Proced Aesthet Dent* 2003; 15: 115-22, 124; quiz 126.
39. **Pieri F, Aldini NN, Fini M, Corinaldesi G.** Immediate occlusal loading of immediately placed implants supporting fixed restorations in completely edentulous arches: A 1-year prospective pilot study. *J Periodontol* 2009; 80: 411-421.
40. **Crespi R, Cappare P, Gherlone E, Romanos GE.** Immediate occlusal loading of implants placed in fresh sockets after tooth extraction. *Int J Oral Maxillofac Implants* 2007; 22: 955-962.
41. **Covani U, Bortolaia C, Barone A, Sbordone L.** Bucco-lingual crestal bone changes after immediate and delayed implant placement. *J Periodontol* 2004; 75: 1605-1612.

42. **Lasella JM, Greenwell H, Miller RL, et al.** Ridge preservation with freeze-dried bone allograft and a collagen membrane compared to extraction alone for implant site development: A clinical and histologic study in humans. *J Periodontol* 2003; 74: 990–999.
43. **Covani U, Cornelini R, Barone A.** Bucco-lingual bone remodeling around implants placed into immediate extraction sockets: A case series. *J Periodontol* 2003; 74: 268–273.
44. **Botticelli D, Berglundh T, Lindhe J.** Hard-tissue alterations following immediate implant placement in extraction sites. *J Clin Periodontol* 2004; 31: 820–828.
45. **Covani U, Cornelini R, Barone A.** Vertical crestal bone changes around implants placed into fresh extraction sockets. *J Periodontol* 2007; 78: 810–815.
46. **Lindeboom JA, Tjiook Y, Kroon FH.** Immediate placement of implants in periapical infected sites: A prospective randomized study in 50 patients. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2006; 101: 705–710.
47. **Crespi R, Cappare P, Gherlone E.** Immediate loading of dental implants placed in periodontally infected and non-infected sites: A 4-year follow-up clinical study. *J Periodontol* 2010; 81: 1140–1146.
48. **Crespi R, Cappare P, Gherlone E.** Fresh-socket implants in periapical infected sites in humans. *J Periodontol* 2010; 81: 378–383.
49. **Pessoa RS, Coelho PG, Muraru L, et al.** Influence of implant design on the biomechanical environment of immediately placed implants: Computed tomography-based nonlinear three-dimensional finite element analysis. *Int J Oral Maxillofac Implants* 2011; 26: 1279–1287.
50. **Peñarrocha-Oltra D, Covani U, Aparicio A, Ata-Ali J, Peñarrocha-Diago M, Peñarrocha-Diago M.** Immediate versus conventional loading for the maxilla with implants placed into fresh and healed extraction sockets to support a full-arch fixed prosthesis restoration: A nonrandomized controlled clinical study. *Int J Oral Maxillofac Implants* 2013; 28: 1116–1124.
51. **Mura P.** Immediate loading of tapered implants placed in postextraction sockets: Retrospective analysis of the 5-year clinical outcome. *Clin Implant Dent Relat Res* 2012; 14: 565–574.
52. **Grunder U, Polizzi G, Goene R, Hatano N, Henry P, Jackson WJ et al.** A 3 year prospective multicenter follow-up report on the immediate and delayed immediate placement of implants. *Int J Oral Maxillofac Implants* 1999; 14: 210–6.
53. **Schwartz-Arad D, Gulayev N, Chashu G.** Immediate versus non-immediate implantation for full-arch fixed reconstruction following extraction of all residual teeth. A retrospective comparative study. *J Periodontol* 2000; 71: 923–8.
54. **Tolman DE, Keller EE.** Endosseous implant placement immediately following dental extraction and alveoloplasty: Preliminary report within 6-year follow-up. *Int J Oral Maxillofac Implants* 1991; 6: 24–8.
55. **Ormianer Z, Piek D, Livne S, et al.** Retrospective clinical evaluation of tapered implants: 10-year follow-up of delayed and immediate placement of maxillary implants. *Implant Dent* 2012; 21: 350–356.
56. **Gelb DA.** Immediate implant surgery: Three-year retrospective evaluation of 50 consecutive cases. *Int J Oral Maxillofac Implants* 1993; 8: 388–399.
57. **Lang NP, Bragger U, Hämmerle CH, Sutter F.** Immediate transmucosal implants using the principle of guided tissue regeneration. I. Rationale, clinical procedures and 30-month results. *Clinical Oral Implants Res* 1994; 5: 154–163.
58. **Vanden Bogaerde L, Rangert B, Wendelhag I.** Immediate/early function of Bränemark System TiUnite implants in fresh extraction sockets in maxillae and posterior mandibles: An 18-month prospective clinical study. *Clin Implant Dent Relat Res* 2005; 7(suppl 1): S121–S130.