LLAVES GALIMPLANT GALIMPLANT WRENCHES CHAVES GALIMPLANT

Español

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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUOVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U. The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations. Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned. Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about its products through the platform: galimplant training.

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

The galimplant® wrenches have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner.

The galimplant® wrenches are reusable. The reusing of these types of products could lead to a deterioration of their characteristics.

The following descriptions are insufficient for the immediate application of any galimplant® wrench, this system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

Description

All wrenches are made with surgical steel and have a 1.31 mm hexagon incorporated. See in the **PRODUCTS** section the galimplant® products compatible with each of the wrenches presented below.

Indications

>Short/long manual wrench

Is used to manually turn and upturn all types of galimplant® screws.

Has a gratory crown in its upper part, which enables the screw to be fixed, and hence reducing the risk of it separating from the wrench.

For added security, this crown has a tiny orifice designed such that the operator can introduce dental floss and hence have complete control over the wrench at all times, avoiding possible aspiration.

>Short/long prosthetic wrench

Is used to manually turn and upturn all types of galimplant® screws.

Permits application of the optimal torque using the torque ratchet, using a direct driver or ratchet adapter.

>Direct short/long prosthetic wrench

Is used to manually turn and upturn all types of galimplant® screws.

Permits application of the optimal torque using the torque ratchet.

> Direct short/long dual prosthetic wrench

Is used to manually turn and upturn all types of galimplant® screws.

Permits application of the optimal torque using the torque ratchet. For added security, this crown has a tiny orifice designed such that the operator can introduce dental floss and hence have complete control over the wrench at all times, avoiding possible aspiration.

>Short/long motorised prosthetic wrench

It's used to manually turn and upturn all types of galimplant® screws.

Is connected directly at the contra-angle enabling the application of adequate torque to the prosthetic screws, supporting a maximum of 25Ncm.

>Motorised wrench

Once connected at the contra-angle, is used to remove the implant from the packaging, transport it too the implant site and introduce it mechanically.

>Manual driver

Is used to manually introduce the implant into its site.

>Direct Short/long driver

Is used to introduce the implant into the site, whilst assisted by the ratchet wrench.

>Direct Overdent wrench

This wrench is designed to be used with retained prosthetics. It has various functions: the square-shaped end, allows manual turning of the overdent abutment onto the implant whilst simultaneous connection of the ratchet wrench permits application of optimal torque (25Ncm). The opposite end fixes the retainer onto the overdent cap (female).

>Direct Short/long multi-position prosthetic wrench

Used to manually turn and upturn all galimplant® straight multi-position abutments.

Allows the application of optimal torque to straight multi-position abutments when using the ratchet wrench.

Patient Groups and Clinical Benefits

The galimplant® original wrench with direct implant connection are intended for use in patients subject to dental implant treatment.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The use of wrenches is not recommended for patients:

- who do not meet the necessary medical requirements for dental implant surgery.
- who are allergic or hypersensitive to medical grade stainless steel.

Secondary effects, interactions and complications

The secondary effects, interactions and complications are those which derive from the implant fitting proceedings.

Cleaning, disinfection and sterilisation

The galimplant® wrenches are sold in a NON-STERILE state. In accordance with the EN ISO 17664 Normative, Nueva Galimplant will provide all of the information necessary to process their medical devices labelled as “Non sterile” and which need to be cleaned, disinfected and sterilized prior to use, both for the first use and for all subsequent usages. All models of the galimplant® wrenches are included in said medical devices.

Warnings:

- The new “non-sterile” products should be processed in the same way as the used products.
- For their own safety, staff should use adequate protective clothing and equipment at all times: gloves, mask, protective glasses, etc;
- Nueva Galimplant recommends the use of aldehyde free cleaner and disinfectants with proven effectiveness, in order to avoid the fixation of blood proteins (included on the VAH/DGHM list or which have the CE Marking).
- The use of neutral or slightly alkaline cleaning agents is recommended by Nueva Galimplant.
- It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.
- Nueva Galimplant cannot control the clinical manipulation procedures, the cleaning methods or any other user level conditions, therefore they will not be held liable for the sterilization of the products even when the recommended guidelines contained in this instruction manual have been followed.

Limitations and restrictions:

The repetition of the indicated processes has a minimal effect on the galimplant® wrenches. The useful life of these is normally determined by the wear and tear and damage caused due to repeated use. The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use.

Preparation in the place of use

As soon as possible following their use, the galimplant® wrenches should be transferred to the place where they are to be cleaned, in order avoid the drying-on of any blood residues, tissues, secretions or bone fragments. If it is likely that the cleaning process will be delayed, these should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination.

Preparation before cleaning

This stage must be carried out regardless as to whether the cleaning and disinfection processes are carried out manually or automatically.

- As soon as possible following their use, immerse the galimplant® wrenches in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Make sure that all the wrench's surfaces are covered by the solution and that the instruments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until all of the visible dirt and residue has been eliminated, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Remove the wrenches from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the wrenches either of the two following methods may be used:

Manual Method

- Place the galimplant® wrenches in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. Nueva Galimplant has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the galimplant® wrenches in the ultrasonic bath and keep them at a temperature of 20° C for 15 minutes.
- Remove the wrenches from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.
- Dry the galimplant® wrenches using filtered compressed air and single use lint free cloth.

Automatic Method

Place the galimplant® wrenches in the thermal washer-disinfector's tray or basket.

Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.

- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Inspection, maintenance and testing

- Carefully inspect all of the galimplant® wrenches in order to check if any damage has been caused to them (corrosion, loss of shape, breaking) or contamination. We recommend that the inspection is carried out in a well-lit environment or using a magnifying lens in order to improve visibility.
- If there are any traces of contamination on any of the devices the cleaning and disinfection process should be repeated.
- In addition, a functional check should be carried out for any devices which have mobile parts in order to check their functionality.

Packaging

Once they are clean, disinfected and completely dry, the wrenches must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.
- Protect the individual devices or surgical boxes sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

Once the devices have been packaged a label should be added to each sterilisation bag which contains data such as the sterilisation date and it should include the reference, batch and the number of devices if the bag does not have a transparent window which allows the user to see what is inside.

Sterilisation

- Load the autoclave

Storage and transport

The galimplant® wrenches are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYMMXXXX).

This packaging contains and protects the product during storage and transport and there is no specified shelf life for it.

The wrenches are made of stable materials under normal environmental conditions, although it is recommended to store them in the original packaging and not expose them to extreme conditions (temperatures, humidity, blows...) that could affect its packaging and labelling.

When removing the product, it must be verified that there are no damages and that the product corresponds to what is described on the label.

WARNING

Galimplant® implants should only be placed using original galimplant® wrenches.

The physician must provide the patient with information on the contraindications, side effects, interactions and complications associated with the galimplant® original wrench.

The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® wrenches and must follow the manufacturer's instructions for use of these materials.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, they should only be used by professionals.

Failure to observe these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Loosening of the abutment or other components.
- Implant failure.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product

The device should be disposed of in accordance with the local laws and the environmental regulations, taking into consideration the different contamination levels.

Additional information

For more information concerning the use of galimplant® products, contact your local representative or visit our website www.galimplant.com.

Validity

These instructions for use will be valid as of January 2022.

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