**LLAVE FIADOR GALIMPLANT®
GALIMPLANT® FIADOR WRENCH
CHAVE FIADOR GALIMPLANT®**

Español

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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical devices. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U. The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations. Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned. Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about it products through the platform: galimplant.training.

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant. The galimplant® fiador wrench has not been sterilised before delivery, so it must be cleaned, disinfected and sterilised before use and it is essential that it is handled in an aseptic manner. The galimplant® fiador wrench is reusable. The reusing of these types of products could lead to a deterioration of its characteristics.

The following descriptions are insufficient for the immediate application of galimplant® fiador wrench. This system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

Description

The galimplant® fiador wrench is made of medical stainless steel. See in the PRODUCTS section the galimplant® products compatible with the wrench presented below.

Indications

- The galimplant® fiador wrench is indicated to perform the following functions:
 - It fixes the implant transfer of the galimplant® implants both to facilitate their removal, and to proceed to modify the depth of insertion of them, occupying the smallest possible occlusal space
 - It allows the friction release of the ratchet wrench adapter, which is provided by the manual screwdriver, to be able to measure both the insertion torque of the implants and to provide the transepithelial abutments and prosthetic screws, with the appropriate torque.
 - The body of the wrench features a series of holes, as well as a pair of raised prominences to distribute and center the implants in fully edentulous upper arches.

Patient Groups and Clinical Benefits

The galimplant® fiador wrench with direct implant connection is intended for use in patients subject to dental implant treatment.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The use of galimplant® fiador wrench is not recommended for patients:

- who do not meet the necessary medical requirements for dental implant surgery.
- who are allergic or hypersensitive to medical grade stainless steel.


Secondary effects, interactions and complications

The secondary effects, interactions and complications are those, which derive from the implant fitting proceedings.

Procedure of use

Step A-B: Center of the prominence in the midline of the maxilla coinciding with the intermaxillary suture.

On both sides of the pair of prominences are two holes 8 mm apart, which allows the correct placement of the implants in positions 11i and 21i. Respecting the biological distances required by both the implants and the intermediate pillars.

Step C: It is rotated on the center of the remaining alveolar rim with respect to the ipsilateral orifice corresponding to the central incisor.

From this first hole (of the central incisor), proceed to mark the next hole, which is at a distance of 12 mm. And, corresponds to the position of the canine (13i or 23i)

Step D: Once the position of the implant bed in positions 13i and / or 23i is marked, the key is moved so that the hole used to locate the central corresponds to that of the canine.

From this second position the key will be rotated on the alveolar process, to go to mark the position of the second premolar (15i, 25i), this distance corresponds to 14 mm, and is located in the second hole that is further away of the one that corresponds with the canine.

Cleaning, disinfection and sterilisation

Attention! The galimplant® fiador wrench is sold in a NON-STERILE state. In accordance with the EN ISO 17664 Normative, Nueva Galimplant will provide all of the information necessary to process their medical devices labelled as "Non sterile" and which need to be cleaned, disinfected and sterilized prior to use, both for the first use and for all subsequent usages. The galimplant® fiador wrench is included in said medical devices.

Warnings:

- The new "non-sterile" products should be processed in the same way as the used products.
- For their own safety, staff should use adequate protective clothing and equipment at all times: gloves, mask, protective glasses, etc;
- Nueva Galimplant recommends the use of aldehyde free cleaner and disinfectants with proven effectiveness, in order to avoid the fixation of blood proteins (included on the VAH/DGHM list or which have the CE Marking).
- Nueva Galimplant recommends the use of neutral or slightly alkaline cleaning agents.
- It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.
- Nueva Galimplant cannot control the clinical manipulation procedures, the cleaning methods or any other user level conditions, therefore they will not be held liable for the sterilization of the products even when the recommended guidelines contained in this instruction manual have been followed.

Limitations and restrictions:

The repetition of the indicated processes has a minimal effect on the galimplant® fiador wrench. The useful life of this one is normally determined by the wear and tear and damage caused due to repeated use. The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use.

Preparation in the place of use

As soon as possible following it use, the galimplant® fiador wrench should be transferred to the place where it is to be cleaned, in order avoid the drying-on of any blood residues, tissues, secretions or bone fragments. If it is likely that the cleaning process will be delayed, this should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination.

Preparation before cleaning

This stage must be carried out regardless as to whether the cleaning and disinfection processes are carried out manually or automatically.

- As soon as possible following its use, immerse the galimplant® fiador wrench in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Makes sure that all the wrench's surface is covered by the solution and that the instruments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until all of the visible dirt and residue has been eliminated, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Removes the wrench from the container containing solution and rinse it thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfects the galimplant® fiador wrench either of the two following methods may be used:

Manual Method

- Places the galimplant® fiador wrench in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. Nueva Galimplant has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the galimplant® fiador wrench in the ultrasonic bath and keep it at a temperature of 20° C for 15 minutes.
- Removes the galimplant® fiador wrench from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.
- Dry the galimplant® fiador wrench using filtered compressed air and single use lint free cloth.

Automatic Method

- Places the galimplant® fiador wrench in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the instrumental we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Inspection, maintenance and testing

- Carefully inspects the galimplant® fiador wrench in order to check if any damage has been caused to it (corrosion, loss of shape, breaking) or contamination. We recommend that the inspection is carried out in a well-lit environment or using a magnifying lens in order to improve visibility.
- If there are any traces of contamination on any of the devices the cleaning and disinfection process should be repeated.
- In addition, a functional check should be carried out for any devices that have mobile parts in order to check their functionality.

Packaging

Once it is clean, disinfected and completely dry, the galimplant® fiador wrench must be packaged in a sterilisation bag. The packaging of thermally sterilised medical devices must comply with the following requisites:

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- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.

- Protect the individual devices or surgical boxes sufficiently from mechanical damage.

- Comply with the EN ISO 11607 normative.

Once the devices have been packaged a label should be added to each sterilisation bag which contains data such as the sterilisation date and it should include the reference, batch and the number of devices if the bag does not have a transparent window which allows the user to see what is inside.

Sterilisation

- Load the autoclave with the conditioned galimplant® fiador wrench in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant recommends using an autoclave, which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.

Storage and transport

The galimplant® fiador wrench is delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYYYMMXXX).

This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

The galimplant® fiador wrench is made of with stable materials under normal environmental conditions, but it is recommended to store it in the original packaging and not expose it to extreme conditions (temperatures, humidity, shock, etc.) that may affect the packaging and labelling.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

WARNING

The galimplant® implants should only be placed using original galimplant® fiador wrench.

The physician must provide the patient with information on the contraindications, side effects, interactions and complications associated with the galimplant® fiador wrench.

The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® fiador wrench and must follow the manufacturer's instructions for use of these materials.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, professionals should only use them.

Failure to observe these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Implant failure.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product

The device should be disposed of in accordance with the local laws and the environmental regulations, taking into consideration the different contamination levels.

Additional information

For more information concerning the use of galimplant® products, contact your local representative or visit our website www.galimplant.com.

Validity

These instructions for use will be valid as of May 2021.

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