PILAR MULTI-POSICIÓN RECTO ESTÉTICO GALIMPLANT
GALIMPLANT AESTHETIC STRAIGHT MULTI-POSITION ABUTMENT
PILAR MULTIPOSIÇÃO RETO ESTÉTICO GALIMPLANT

Español

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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned. Nueva Galimplant with the aim of ensuring the correct use of its products offers courses and practical sessions about it products through the platform: galimplant.training.

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

Galimplant® aesthetic straight multi-position abutments have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner.

Galimplant® aesthetic straight multi-position abutments are of single use. The reusing of these types of products could lead to a deterioration of their characteristics which involves risk of tissue infection, surgical or prosthodontic failure and/or damage to the patient's health.

The following descriptions are insufficient for the immediate application of any abutment from the galimplant™ system, this system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

Description

Galimplant® aesthetic straight multi-position abutments are made from grade 5 titanium and can be used in galimplant® implants with different connection types (external and internal), platforms, diameters and lengths.

See in the PRODUCTS section the implant models compatible with each of the configurations presented below.

Our straight aesthetic abutments are available in the following configurations:

- Internal Connection
 - Rotational aesthetic straight multi-position abutment with shoulder heights of 1, 2, 3, 4 and 5 mm for narrow platform and heights of 2, 3, 4 and 5mm for wide platform and pink colour anodized.
 - Non-rotational aesthetic straight multi-position abutment with shoulder heights of 1, 2, 3, 4 and 5 mm for narrow platform and golden colour anodized.
 - Solidary non-rotational aesthetic straight multi-position abutment with shoulder heights of 1, 2, 3, 4 and 5 mm for narrow platform and heights of 2, 3, 4 and 5mm for wide platform and golden colour anodized.
 - Rotational (pink colour anodized) or non-rotational (golden colour anodized) slim straight aesthetic multi-position abutment with shoulder heights of 2, 3, 4, 5 and 6 mm.
 - Transformable to external aesthetic straight multi-position abutment (golden anodized) with shoulder height 2, 3 and 4mm
- External Connection
 - Non-rotational aesthetic straight multi-position abutment with shoulder heights of 2, 3 and 5 mm and golden colour anodized.
 - May be converted to rotational using the octagonal castable abutment.

The shoulder height will depend on the soft tissue height.

They come with:

-2 screws blue colour anodized (one clinical, packed inside a capsule and another prosthetic).

-1 transepithelial transfer attachment

-1 screw for the transepithelial transfer attachment

-1 transepithelial closing abutment

Recommended torque:

-30 Ncm aesthetic straight multi-position abutment.

-25 Ncm prosthetic screw.

Do not reuse. Single use.

Intended use

Our aesthetic straight multi-position abutments with direct implant connection are designed for screw-retained restorations, with the aim of restoring the patient's mastication function and phonetics.

Indications

Our aesthetic straight multi-position abutments are recommended in fixed single- or multi-unit screw-retained restorations with high aesthetic requirements.

Our non-rotational multi-position slim straight aesthetic abutments are not recommended for large fixed multi-unit screw-retained restorations.

Patient Groups and Clinical Benefits

The galimplant® aesthetic straight multi-position abutments with direct implant connection are intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that require treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use. The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The use of aesthetic straight multi-position abutments is not recommended for patients:

- who do not meet the necessary medical requirements for dental implant surgery.
- who are allergic or hypersensitive to commercially pure titanium or Ti-6Al-4V titanium alloy (titanium, aluminium and vanadium).

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures.

Possible complications following prosthetic restoration include:

- Complications: Loosening/fracture of the abutment and / or screw, mechanical problems associated with the prosthesis, loss of prosthesis retention, loss of the hole sealing of the screw-retained prosthesis and microgap in the implant-abutment connection.

If one of these disorders/complications appears, the patient should consult their specialist doctor who will give them the corresponding instructions.

Storage and transport

The aesthetic straight multi-position abutments are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYYYMMXXX). This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

The aesthetic straight multi-position abutments must be stored in a dry place in the original packaging at room temperature and not exposed to direct sunlight.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

Cleaning/Sterilisation/Conservation

Caution! All our aesthetic straight multi-position abutments are sold in a NON-STERILE state. Before use, they must be cleaned, disinfected and sterilised, observing the following procedures validated by Nueva Galimplant SLU. Said procedures should be carried out before intraoral use.

The repetition of the procedures described in this section do not alter the characteristics of these devices. The non-observance of these instructions may cause cross infections.

The use of neutral or slightly alkaline cleaning agents is recommended. It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.

Preparation before cleaning

- Immerse all the attachments in a container that contains a cleaning solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time.

- Make sure that all the attachments surface are covered by the solution and that the attachments do not touch each other.

- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until eliminating all the visible dirt and residue, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.

- Remove the attachments from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the attachments either of the two following methods may be used:

Manual Method

- Place the attachments in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. *Nueva Galimplant S.L.U. has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.*
- Immerse the basket with the attachments in the ultrasonic bath and keep them at a temperature of 20° C.
- Remove the attachments from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.

Automatic Method

- Place the attachments in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. *Nueva Galimplant S.L.U. has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:*
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Packaging

NUEVA GALIMPLANT,S.L.U

C/ Benigno Quiroga, 90 27600 SARRIA LUGO ESPAÑA