**INTERFASES PARA PILARES MULTI-POSICIÓN GALIMPLANT
GALIMPLANT INTERFACES FOR MULTI-POSITION ABUTMENTS
INTERFACES PARA PILARES MULTIPOSIÇÃO GALIMPLANT**

Español

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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

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application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U. The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to

Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.

- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant S.L.U. has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55°C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93°C thermal disinfection for 5 minutes.
 - 110°C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfecter manufacturers.

Packaging

Once they are clean, disinfected and completely dry, the attachments must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137°C and be sufficiently permeable for steam.
- Protect the devices sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

Sterilization

Load the autoclave with the conditioned attachments in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.

Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.

Nueva Galimplant S.L.U. recommends using an autoclave which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance.

The validation of the sterilisation was carried out in a steam autoclave using a 134°C pre-vacuum cycle for 6 minutes.

To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.

Caution! The re-use and/or re-sterilisation of prosthetic abutments can cause risk of infection and/or cross contamination.

Documentation, traceability and implant card.

Nueva Galimplant S.L.U. recommends keeping an exhaustive clinical, radiological and photographic documentation. All products can be identified through their reference and batch code and / or UDI code. To maintain the mandatory traceability of these products, two adhesive labels are included, one for the medical record and the other for the patient, the latter must be attached to the implant card supplied with the products and which must be delivered to the patient fully completed once the surgery is over.

WARNINGS

The galimplant® interfaces for multi-position abutments should be placed only on multi-positions straight and angled original abutments.

The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® interfaces for multi-position abutments and must follow the manufacturer's instructions for use of these materials.

Before the prosthetic part, the physician should carry out an exhaustive patient study, in order to determine a correct planning of the case.

For the placement of galimplant® interfaces for multi-position abutments, use only original wrenches of the brand galimplant®. Consult the compatible product models in the catalogue available at www.galimplant.com.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, they should only be used by professionals.

The prosthetic abutments must solely and exclusively be placed during the implants osteointegration or if the implants have a great primary stability.

Failure to observe the procedures set out in these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Implant failure.

To guarantee a good long term treatment result, the clinician must inform the patient about appropriate home oral hygiene and the importance of attending periodic clinical and radiographic consultations to perform a control and preventive diagnosis of any side effect and/or complication.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

The practitioner must provide the patient with the information about side effects, interactions and complications in relation galimplant® interfaces for multi-position abutments.

Procedures for the disposal of the product

Once they have been extracted from the oral cavity, the prosthetic components must be disposed of as biological waste, in accordance with current local laws.

Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) for the product is available in the European database EUDAMED. The document can be obtained from the website: <https://ec.europa.eu/tools/eudamed>.

Validity

These instructions for use will be valid as of May 2022.

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