SISTEMA DE REGISTRO OCLUSAL GALIMPLANT®

GALIMPLANT® OCCLUSAL REGISTRATION SYSTEM

SISTEMA DE REGISTRO OCLUSAL GALIMPLANT®

Español

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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical devices. The Galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the Galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with Galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about it products through the platform: galimplant.training.

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

The galimplant® occlusal registration system is delivered without sterilization, so it is necessary to clean, disinfect and sterilize it before use and aseptic aseptic handling is essential.

The galimplant® occlusal registration system is reusable. The reuse of this type of products leads to a possible deterioration of its characteristics.

The following descriptions are not sufficient for the immediate use of the galimplant® occlusal registration system, this system should only be used by dentists, doctors and surgeons. The dentist / physician / surgeon must have knowledge in oral rehabilitation and experience in the manipulation of the galimplant® product, to be able to make a safe and correct use of it, following the indications of these instructions of use.

Description

All the pieces that make up the occlusal registration system are made of surgical steel.

This system is composed of 6 pieces:

- Screwdriver
- Intermaxillary fixator
- Central pointer (vertical dimension regulator)
- Mandibular plate horizontal displacement (large)
- Mandibular plate horizontal displacement (small)
- Rule to determine the vertical dimension

Indications

The occlusal registration system is designed to be used to determine the relationship between the jaws.

Patient Groups and Clinical Benefits

The galimplant® occlusal registration system with direct implant connection are intended for use in patients subject to dental implant treatment.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The use of the occlusal registration system is not recommended for patients:

- who are allergic or hypersensitive to surgical steel.

Clinical Steps

The upper and lower acrylic plates are positioned on the alveolar rims. The central pointer is adjusted to determine the vertical dimension (the millimeter rule is used to measure the golden proportions, taking as reference nasion, subnasal, pogonion). The patient is instructed to close the mouth and contact the central pointer with the horizontal mandibular plate. The patient will perform protrusive movements, right and left laterality, repeatedly and freely (without manipulation). Then remove the bottom plate, mark the surface of the plate with an indelible marker and repeat the movements. The ball will draw the gothic arch on the pen ink in the shape of an arrowhead. Remove the lower plate and fix the intermaxillary fixator with the rear screw of the plate, making the point of the arrow marked on the lower plate coincide with the center of the sphere anchor. Fix the upper and lower plates by fitting the ball into the hole of the intermaxillary fixator. Firmly fix the front screw of the intermaxillary fixator with the screwdriver. Mount in articulator.

Cleaning, disinfection and sterilization

The galimplant® occlusal registration system is sold in a NON-STERILE state. In accordance with the EN ISO 17664 Normative, Nueva Galimplant will provide all of the information necessary to process their medical devices labelled as "Non sterile" and which need to be cleaned, disinfected and sterilized prior to use, both for the first use and for all subsequent usages. The galimplant® occlusal registration system is included in said medical devices.

IMPORTANT! The box containing the pieces of the galimplant® occlusal registration system is not suitable for cleaning, disinfection and sterilization processes, so both processes must be carried out outside the box.

Warnings:

- The new "non-sterile" products should be processed in the same way as the used products.
- For their own safety, staff should use adequate protective clothing and equipment at all times: gloves, mask, protective glasses, etc;
- Nueva Galimplant recommends the use of aldehyde free cleaner and disinfectants with proven effectiveness, in order to avoid the fixation of blood proteins (included on the VAH/DGHH list or which have the CE Marking).
- Nueva Galimplant recommends the use of neutral or slightly alkaline cleaning agents.
- It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.
- Nueva Galimplant cannot control the clinical manipulation procedures, the cleaning methods or any other user level conditions, therefore they will not be held liable for the sterilization of the products even when the recommended guidelines contained in this instruction manual have been followed

Limitations and restrictions

The repetition of the indicated processes has a minimal effect on the galimplant® occlusal registration system. The useful life of these is normally determined by the wear and tear and damage caused due to repeated use.

The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use.

Preparation at the place of use

As soon as possible following their use, the galimplant® occlusal registration system should be transferred to the place where they are to be cleaned, in order avoid the drying-on of any blood residues or secretions. If it is likely that the cleaning process will be delayed, these should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination.

Preparation before cleaning

This stage must be carried out regardless as to whether the cleaning and disinfection processes are carried out manually or automatically.

- As soon as possible following their use, immerse the pieces of the galimplant® occlusal registration system in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Make sure that all the pieces's surfaces are covered by the solution and that the instruments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until all of the visible dirt and residue has been eliminated, paying special attention to the zones that may not be reached by the brushing action, such as grooves or internal parts.
- Remove the pieces from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the pieces of the galimplant® occlusal registration system either of the two following methods may be used:

Manual Method

- Place the pieces in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. Nueva Galimplant has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the pieces in the ultrasonic bath and keep them at a temperature of 20° C for 15 minutes.
- Remove the pieces from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.
- Dry the pieces using filtered compressed air and single use lint free cloth.

NUEVA GALIMPLANT,S.L.U

C/ Benigno Quiroga, 90 27600 SARRIA LUGO ESPAÑA

Automatic Method

- Place the pieces in the thermal washer-disinfector’s tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the instrumental we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Inspection, maintenance and testing

- Carefully inspect all of the Galimplant® occlusal registration system in order to check if any damage has been caused to them (corrosion, loss of shape, breaking) or contamination. We recommend that the inspection is carried out in a well-lit environment or using a magnifying lens in order to improve visibility.
- If there are any traces of contamination on any of the devices the cleaning and disinfection process should be repeated.
- In addition, a functional check should be carried out for any devices that have mobile parts in order to check their functionality.

Packing

Once they are clean, disinfected and completely dry, the pieces must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.
- Protect the individual devices or surgical boxes sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

Once the devices have been packaged a label should be added to each sterilisation bag that contains data such as the sterilisation date and it should include the reference, batch and the number of devices if the bag does not have a transparent window that allows the user to see what is inside.

Sterilization

• Load the parts that make up the galimplant® occlusal registration system conditioned in the sterilization bags in the autoclave. Do not overload the autoclave and always follow the equipment manufacturer’s instructions.

Note: the box containing the parts of the occlusal registration system is not sterilizable, so the parts must be packaged and sterilized out of the box.

- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant recommends using an autoclave that complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance.

The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer’s instructions and keep maintenance up to date.

Storage and transport

The galimplant® occlusal registration system is delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYYYMMXXX). This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

The parts that make up the occlusal registration system are made of with stable materials under normal environmental conditions, but it is recommended to store them in the original packaging and not expose them to extreme conditions (temperatures, humidity, shock, etc.) that may affect the packaging and labeling.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

WARNINGS

The physician must provide the patient with information on the contraindications, side effects, interactions and complications associated with the galimplant® occlusal registration system.

The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® occlusal registration system and must follow the manufacturer's instructions for use of these materials.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, professionals should only use them.

Failure to observe these instructions could harm the patient and give rise to one or all of the following complications:

- Infection.
- Damage to the components or instruments.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product

The device should be disposed of in accordance with the local laws and the environmental regulations, taking into consideration the different contamination levels.

Additional information

For more information concerning the use of galimplant® products, contact your local representative or visit our website www.galimplant.com.

Validity

These instructions for use will be valid as of May 2021.

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