IMPLANTE DENTAL ESPECIAL ONE-PIECE
ONE-PIECE SPECIAL DENTAL IMPLANT
IMPLANTE DENTÁRIO ESPECIAL ONE-PIECE

Español

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English

IMPORTANT! Read the information included in these instructions for use before using the product

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The Galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the Galimplant® System is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions.

This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their

obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about it products through the platform: galimplant training.

Precautions and warnings

The galimplant® dental implants are delivered sterile; they are for single use and must not be used after their expiration date.

Warning! Re-sterilisation and/or reuse of implants may result in the risk of crossed contamination.

The following descriptions are insufficient for the immediate application of any Galimplant® system implants. This system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of the product in order to use them safely and correctly following the indications of these instructions for use.

Description

The one-piece implants are manufactured using commercially pure grade IV titanium, in accordance with directive ISO 5832-2 and are mechanized in one piece (implant body, collar and retention mechanism or abutment). The threaded body has a Nanoblastplus® surface, the result of a technique which combines a jet of sand and a triple acid etching, the collar is concave and polished and at coronal part they have a machining retention mechanism with a geometrical figure.

Provided with a titanium implant transfer fixed to the retention mechanism. The implant transfer is manufactured using grade V titanium (Ti6Al4V) in accordance with ISO 5832-3.

The one-piece special implants are available with two types of configurations:

- **Mini implant for cementing (IMC):** the retention mechanism at the coronal part has a hexagonal figure.
- **Overdenture implant (ISD):** the retention mechanism at the coronal part has an overdent type geometric figure.

The selection of the model, diameter and length will depend on the height and width of the bone, the bone density and the surgical technique which is to be used. Check out the implant models available in the PRODUCTS section.

Compatibility

The IMC implant model is designed for use with the IMC castable abutment. The ISD implant model is compatible with the overdent system.

Galimplant® surgical drill system must be used to prepare the implant bed, following the indicated drilling system and only the original galimplant® wrenches should be used for the handling and placement of implants, as stated in the instructions for use.

Intended use

The galimplant® implant system is designed for use in oral implantology in order to re-establish lost functions.

Indications

- **Mini implant for cementing (IMC):** are indicated for fixed cemented prosthesis at compact spaces.
- **Overdenture implant (ISD):** are indicated to retain/support removable prosthesis.

Patient Groups and Clinical Benefits

The galimplant® implant system is intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The applicable contraindications are those related to surgical procedures for the fitting of implants.

-Absolute contraindications: Severe internal medical problems, metabolic bone disorders, uncontrolled coagulopathy, un-cooperative or un-motivated patients, alcoholism or drug addiction, psychosis, history of treatment resistant disorder, xerostomia, weakened immune system, illnesses treated with regular usage of steroids, allergy or hypersensitivity to titanium, uncontrollable endocrine illnesses.

-Relative contraindications, in these specific circumstances in which the risk may not be considered as high as the benefit obtained after implant placement, the clinician must assess these conditions: Previously irradiated bone, diabetes, anticoagulation therapy/hemorrhagic diathesis, bruxism, parafunctional habits, unfavourable bone anatomy, smoking, uncontrolled periodontitis, temporomandibular joint disorder, maxillary disease and altered oral mucous membrane susceptible to treatment, controllable endocrine illnesses, pregnancy, lactation, children with unfinished bone growth, insufficient oral hygiene, insufficient bone available or inadequate bone quality, localized radicular remains.

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures. Possible complications following dental implant placement procedures include:

- Temporary disorders: Pain, inflammation, difficulty speaking, gingivitis, edema, swelling, difficulty opening the mouth, bad breath, labial laceration, haemorrhage, contusion, hematoma, temporal paraesthesia.
- Longer- lasting problems: Chronic pain associated with the dental implant, permanent paraesthesia, dyesthesia, bone loss from the maxillary crest, localised or systemic infection, oral-antral or oral-nasal fistulas, periimplantitis, implant failure, irreversible damage to neighbouring teeth, implant, maxillary or prosthetic fractures, aesthetic problems, nerve damage, exfoliation, hyperplasia, periapical infection, mucositis, marginal bone loss, sinusitis, perforation of the anatomical areas.

In case in which any of these disorders appear, the patient should report them to the specialist who will give the corresponding indications.

Packaging and sterility

The one-piece special dental implants have been manufactured and packaged in a controlled environment prior to their sterilisation by gamma rays.

The metallic implant transfer that holds the endosseous implant is presented assembled on a plastic support. The aforementioned plastic support is presented inside a container that guarantees sterility by means of an operculum.

Both the external and the internal packaging bears a label which indicates the following concepts: implant model, size, batch and expiration date, in addition the external label seals the packaging guaranteeing that the content cannot be manipulated.

Nueva Galimplant S.L.U. decline all responsibility in the case of re-sterilizing implants, irrespective of the person performing the sterilization or the method used. A previously used or non-sterile implant should not be implanted under any circumstance. Nueva Galimplant S.L.U. cannot accept the return of an implant if the packaging is damaged.

Storage and transport

The implants should be maintained in their original packaging, in a dry environment, at room temperature and away from direct sunlight. These same conditions are also required during transportation and in addition the product should not be knocked or crushed.

The intact sterile packaging protects the implant from any external influences during storage and transportation, and if it is stored in adequate conditions it guarantees the sterility until the expiration date. The implant should not be used following the expiration date (printed on the label).

When taking the product out, ensure that there is no damage and that the product is as described on the label.

Planning procedures for use

Planning

Before surgery, an exhaustive patient evaluation and study should be carried out (assessment, study models, panoramic and periapical x-rays and C.A.T. scan) in order to determine the anatomical characteristics, occlusion conditions, periodontal condition and to determine the bone suitability in order to draw up the treatment plan. The physician is exclusively authorised to plan this treatment. Faults in the treatment plan may cause the loss of implants.

Procedure

- **Preparation of the implant bed**

The implant bed should be treated in a non-aggressive manner in order to create optimum conditions for the cicatrization of the implant. A heat trauma will hinder or prevent the healing of the dental implant. For this reason, the following measures should be applied in order to reduce, as far as possible, the prevalence of excessive temperatures:

- Use the appropriate drilling speed for each drill diameter. A speed of between 300 and 800 rpm is recommended.
- Use sharpened surgical instruments.
- An intermittent drilling technique should be applied.
- The drills should be highly refrigerated using pre-cooled (5°) sterile physiological saline solution (NaCl).
- Use drills with increasing diameter: first the lance drill and then the Ø 2mm drill. To favour the primary stability of the implant, the diameter of the last drill should be less than the diameter of the implant.

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• **Extraction of the implant from the packaging and insertion into the bed.**
The placement of the implant once the implant bed has been prepared may be done manually, using the manual wrench or mechanically using the contra-angle or the torque wrench. When mechanically placing one-piece special implants, the speed should not exceed 30 rpm and the insertion torque should be between 25-30Ncm.
To take the implant out of the container open the external packaging, unscrew the plastic top and remove the operculum. The implant may be taken out of the package in the following ways:

1. Connect the manual driver to the metal implant transfer and extract the plastic holder from the container.
2. Insert the motor wrench in the contra-angle, connect the point to the metal implant transfer and extract the plastic holder from the container.

Once connected remove the implant from the plastic holder by applying light pressure in an upwards direction at the aperture where it is placed and hold it with any of the assembly techniques indicated above to bring it to the mouth and insert it into the implant bed. Once the implant has been inserted, the implant transfer must be removed by pulling gently.

Documentation, traceability and implant card.
Nueva Galimplant S.L.U. recommends keeping an exhaustive clinical, radiological and photographic documentation. All products can be identified through their reference and batch code and / or UDI code. To maintain the mandatory traceability of these products, two adhesive labels are included, one for the medical record and the other for the patient, the latter must be attached to the implant card supplied with the products and which must be delivered to the patient fully completed once the surgery is over.

WARNINGS
The physician should provide the patient with information regarding the contraindications, side effects, interactions and complications associated with the galimplant® implants.
The patient should avoid activities which require a high level of physical effort immediately after the insertion of the dental implant.
When using our products, the necessary measures should be taken to avoid inhalation. Before placing the implant ensure that there is no damage has been caused to both containers. A damaged container may mean that the implant is now not sterile. When removing the implant from the sterile container follow the corresponding aseptic principles. The container should only be opened when the surgical intervention is to be carried out.
The health professional must ensure that the materials and techniques that are to be used during the clinical and prosthetic procedures are appropriate and compatible with the one-piece special implants and the instructions for use provided by the manufacturer of these materials should be followed.
For the preparation of the implant bed, and the placement and restoration of the one-piece special implants only original galimplant® products must be used.
Failure to follow the procedures summarised in these instructions may cause damage to the patient and may to one or more of the following complications.

- Inhalation or swallowing of the components.
- Infection
- Damage to the implant, abutment, components or instruments.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product
Once extracted from the oral cavity the prosthetic components must be disposed of as biological waste in accordance with the local regulations currently in force.

Summary of Safety and Clinical Performance
The Summary of Safety and Clinical Performance (SSCP) for the product is available in the European database EUDAMED. The document can be obtained from the website: <https://ec.europa.eu/tools/eudamed>.

Validity
These instructions for use will be valid as of May 2022.

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